

Product datasheet for **GA106118**

Human AP4M1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	AP4M1 CRISPRa kit - CRISPR gene activation of human adaptor related protein complex 4 subunit mu 1
Symbol:	AP4M1
UniProt ID:	<u>O00189</u>
Synonyms:	CPSQ3; MU-4; MU-ARP2; SPG50
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)