

## Product datasheet for **GA106130**

### Human **KCNAB3** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	KCNAB3 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily A regulatory beta subunit 3
Symbol:	KCNAB3
UniProt ID:	<u><a href="#">O43448</a></u>
Synonyms:	AKR6A9; KCNA3.1B; KCNA3B; KV-BETA-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)