

Product datasheet for **GA106143**

Human XPR1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	XPR1 CRISPRa kit - CRISPR gene activation of human xenotropic and polytropic retrovirus receptor 1
Symbol:	XPR1
UniProt ID:	<u>Q9UBH6</u>
Synonyms:	SYG1; X3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)