

Product datasheet for **GA106151**

Human **MAGI1** activation kit by **CRISPRa**

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | MAGI1 CRISPRa kit - CRISPR gene activation of human membrane associated guanylate kinase, WW and PDZ domain containing 1 |
| Symbol: | MAGI1 |
| UniProt ID: | Q96QZ7 |
| Synonyms: | AIP-3; AIP3; BAIAP1; BAP-1; BAP1; MAGI-1; Magi1d; TNRC19; WWP3 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)