

Product datasheet for **GA106176**

Human PDLIM7 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PDLIM7 CRISPRa kit - CRISPR gene activation of human PDZ and LIM domain 7
Symbol:	PDLIM7
Components:	GA106176G1, PDLIM7 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106176G2, PDLIM7 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106176G3, PDLIM7 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9NR12
Synonyms:	LMP1; LMP3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)