

Product datasheet for **GA106183**

Human ITGB1BP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ITGB1BP1 CRISPRa kit - CRISPR gene activation of human integrin subunit beta 1 binding protein 1
Symbol:	ITGB1BP1
Components:	GA106183G1, ITGB1BP1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106183G2, ITGB1BP1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106183G3, ITGB1BP1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q14713</u>
Synonyms:	ICAP-1A; ICAP-1alpha; ICAP-1B; ICAP1; ICAP1A; ICAP1B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)