

Product datasheet for **GA106191**

Human GPR37L1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	GPR37L1 CRISPRa kit - CRISPR gene activation of human G protein-coupled receptor 37 like 1
Symbol:	GPR37L1
Components:	GA106191G1, GPR37L1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106191G2, GPR37L1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106191G3, GPR37L1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O60883</u>
Synonyms:	ET(B)R-LP-2; ETBR-LP-2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)