

Product datasheet for **GA106203**

Human ASIC3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ASIC3 CRISPRa kit - CRISPR gene activation of human acid sensing ion channel subunit 3
Symbol:	ASIC3
UniProt ID:	<u>Q9UHC3</u>
Synonyms:	ACCN3; DRASIC; SLNAC1; TNaC1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)