

## Product datasheet for **GA106204**

### Human Kv2.2 (KCNB2) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	KCNB2 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel subfamily B member 2
Symbol:	KCNB2
UniProt ID:	<a href="#">Q92953</a>
Synonyms:	KV2.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)