

Product datasheet for **GA106207**

Human P311 (NREP) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	NREP CRISPRa kit - CRISPR gene activation of human neuronal regeneration related protein
Symbol:	NREP
UniProt ID:	Q16612
Synonyms:	C5orf13; D4S114; P311; PRO1873; PTZ17; SEZ17
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)