

Product datasheet for **GA106230**

Human TAOK2 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | TAOK2 CRISPRa kit - CRISPR gene activation of human TAO kinase 2 |
| Symbol: | TAOK2 |
| UniProt ID: | Q9UL54 |
| Synonyms: | MAP3K17; PSK; PSK1; PSK1-BETA; TAO1; TAO2 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)