

## Product datasheet for **GA106241**

### Human PPIG activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | PPIG CRISPRa kit - CRISPR gene activation of human peptidylprolyl isomerase G |
| Symbol:       | PPIG  |
| UniProt ID:   | <a href="#">Q13427</a>  |
| Synonyms:     | CARS-Cyp; CYP; SCAF10; SRCyp  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)       |

#### Product images:



[View online >](#)