

Product datasheet for **GA106253**

Human PLAA activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLAA CRISPRa kit - CRISPR gene activation of human phospholipase A2 activating protein
Symbol:	PLAA
Components:	GA106253G1, PLAA gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106253G2, PLAA gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106253G3, PLAA gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y263</u>
Synonyms:	DOA1; PLA2P; PLAP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)