

## Product datasheet for **GA106280**

### Human SNRNP40 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SNRNP40 CRISPRa kit - CRISPR gene activation of human small nuclear ribonucleoprotein U5 subunit 40
Symbol:	SNRNP40
UniProt ID:	<u><a href="#">Q96DI7</a></u>
Synonyms:	40K; HPRP8BP; PRP8BP; PRPF8BP; SPF38; WDR57
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)