

Product datasheet for **GA106284**

Human TJP2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TJP2 CRISPRa kit - CRISPR gene activation of human tight junction protein 2
Symbol:	TJP2
UniProt ID:	Q9UDY2
Synonyms:	C9DUPq21.11; DFNA51; DUP9q21.11; PFIC4; X104; ZO2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)