

Product datasheet for **GA106292**

Human KCNK6 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNK6 CRISPRa kit - CRISPR gene activation of human potassium two pore domain channel subfamily K member 6
Symbol:	KCNK6
UniProt ID:	Q9Y257
Synonyms:	K2p6.1; KCNK8; TOSS; TWIK-2; TWIK2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)