

Product datasheet for **GA106305**

Human QK1 (QKI) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	QKI CRISPRa kit - CRISPR gene activation of human QKI, KH domain containing RNA binding
Symbol:	QKI
UniProt ID:	<u>Q96PU8</u>
Synonyms:	Hqk; hqkl; QK; QK1; QK3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)