

Product datasheet for **GA106340**

Human PIGL activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIGL CRISPRa kit - CRISPR gene activation of human phosphatidylinositol glycan anchor biosynthesis class L
Symbol:	PIGL
UniProt ID:	Q9Y2B2
Synonyms:	CHIME
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)