

Product datasheet for **GA106341**

Human PIGB activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIGB CRISPRa kit - CRISPR gene activation of human phosphatidylinositol glycan anchor biosynthesis class B
Symbol:	PIGB
UniProt ID:	Q92521
Synonyms:	GPI-MT-III; PIG-B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)