

Product datasheet for GA106343

Human PSMF1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMF1 CRISPRa kit - CRISPR gene activation of human proteasome inhibitor subunit 1
Symbol:	PSMF1
UniProt ID:	<u>Q92530</u>
Synonyms:	PI31
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)