

Product datasheet for **GA106362**

Human LITAF activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LITAF CRISPRa kit - CRISPR gene activation of human lipopolysaccharide induced TNF factor
Symbol:	LITAF
UniProt ID:	Q99732
Synonyms:	PIG7; SIMPLE; TP53I7
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)