

Product datasheet for **GA106384**

Human PIG3 (TP53I3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TP53I3 CRISPRa kit - CRISPR gene activation of human tumor protein p53 inducible protein 3
Symbol:	TP53I3
UniProt ID:	Q53FA7
Synonyms:	PIG3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)