

## Product datasheet for **GA106387**

### Human IGDCC3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	IGDCC3 CRISPRa kit - CRISPR gene activation of human immunoglobulin superfamily DCC subclass member 3
Symbol:	IGDCC3
UniProt ID:	<a href="#">Q8IVU1</a>
Synonyms:	HsT18880; PUNC
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)