

Product datasheet for **GA106423**

Human WTAP activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | WTAP CRISPRa kit - CRISPR gene activation of human WT1 associated protein |
| Symbol: | WTAP |
| UniProt ID: | Q15007 |
| Synonyms: | Mum2 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)