

Product datasheet for **GA106443**

Human Kallikrein 4 (KLK4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KLK4 CRISPRa kit - CRISPR gene activation of human kallikrein related peptidase 4
Symbol:	KLK4
Components:	GA106443G1, Kallikrein 4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106443G2, Kallikrein 4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106443G3, Kallikrein 4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y5K2</u>
Synonyms:	AI2A1; ARM1; EMSP; EMSP1; kallikrein; KLK-L1; PRSS17; PSTS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)