


## Product datasheet for **GA106444**

### Human **TCL1B** activation kit by **CRISPRa**

#### Product data:

|               |                                                                                                                                                                                                                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| Description:  | TCL1B CRISPRa kit - CRISPR gene activation of human T cell leukemia/lymphoma 1B                                                                                                                                                                                                                            |
| Symbol:       | TCL1B                                                                                                                                                                                                                                                                                                      |
| Components:   | GA106444G1, TCL1B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA106444G2, TCL1B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA106444G3, TCL1B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">O95988</a>                                                                                                                                                                                                                                                                                     |
| Synonyms:     | SYN-1; TML1                                                                                                                                                                                                                                                                                                |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:



[View online »](#)