

Product datasheet for **GA106461**

Human FISH (SH3PXD2A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SH3PXD2A CRISPRa kit - CRISPR gene activation of human SH3 and PX domains 2A
Symbol:	SH3PXD2A
UniProt ID:	Q5TCZ1
Synonyms:	FISH; SH3MD1; TKS5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)