

Product datasheet for **GA106468**

Human PLCH2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLCH2 CRISPRa kit - CRISPR gene activation of human phospholipase C eta 2
Symbol:	PLCH2
Components:	GA106468G1, PLCH2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106468G2, PLCH2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106468G3, PLCH2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O75038</u>
Synonyms:	PLC-eta2; PLC-L4; PLCeta2; PLCL4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)