

## Product datasheet for **GA106469**

### Human **TTC37** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	TTC37 CRISPRa kit - CRISPR gene activation of human tetratricopeptide repeat domain 37
Symbol:	TTC37
UniProt ID:	<u><a href="#">Q6PGP7</a></u>
Synonyms:	KIAA0372; Ski3; THES
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)