

Product datasheet for **GA106474**

Human IQ motif containing B1 (IQCB1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IQCB1 CRISPRa kit - CRISPR gene activation of human IQ motif containing B1
Symbol:	IQCB1
Components:	GA106474G1, IQ motif containing B1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106474G2, IQ motif containing B1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106474G3, IQ motif containing B1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q15051
Synonyms:	NPHP5; PIQ; SLSN5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)