

Product datasheet for **GA106479**

Human LIMKAIN B1 (KIAA0430) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MARF1 CRISPRa kit - CRISPR gene activation of human meiosis regulator and mRNA stability factor 1
Symbol:	MARF1
Components:	GA106479G1, LIMKAIN B1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106479G2, LIMKAIN B1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106479G3, LIMKAIN B1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9Y4F3
Synonyms:	LKAP; MARF1; PPP1R34
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)