

Product datasheet for **GA106483**

Human EIF5B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EIF5B CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 5B
Symbol:	EIF5B
Components:	GA106483G1, EIF5B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106483G2, EIF5B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106483G3, EIF5B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O60841</u>
Synonyms:	IF2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)