


Product datasheet for **GA106486**

Human Syndecan 3 (SDC3) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | SDC3 CRISPRa kit - CRISPR gene activation of human syndecan 3 |
| Symbol: | SDC3 |
| UniProt ID: | <u>O75056</u> |
| Synonyms: | SDCN; SYND3 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)