

Product datasheet for **GA106490**

Human HISPPD2A (PPIP5K1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PPIP5K1 CRISPRa kit - CRISPR gene activation of human diphosphoinositol pentakisphosphate kinase 1
Symbol:	PPIP5K1
Components:	GA106490G1, HISPPD2A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106490G2, HISPPD2A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106490G3, HISPPD2A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6PFW1</u>
Synonyms:	DKFZp313L0221; HISPPD2A; IP6K; IPS1; KIAA0377; MGC51871; VIP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)