

Product datasheet for **GA106515**

Human DHX34 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DHX34 CRISPRa kit - CRISPR gene activation of human DExH-box helicase 34
Symbol:	DHX34
Components:	GA106515G1, DHX34 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106515G2, DHX34 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106515G3, DHX34 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q14147</u>
Synonyms:	DDX34; HRH1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)