

Product datasheet for **GA106535**

Human RAB11FIP3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RAB11FIP3 CRISPRa kit - CRISPR gene activation of human RAB11 family interacting protein 3
Symbol:	RAB11FIP3
UniProt ID:	<u>Q75154</u>
Synonyms:	CART1; Rab11-FIP3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)