

Product datasheet for **GA106537**

Human KIAA0408 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KIAA0408 CRISPRa kit – CRISPR gene activation of human KIAA0408
Symbol:	KIAA0408
Components:	GA106537G1, KIAA0408 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106537G2, KIAA0408 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106537G3, KIAA0408 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZU52</u>
Synonyms:	RP3-403A15.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

