

Product datasheet for **GA106541**

Human SART3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SART3 CRISPRa kit - CRISPR gene activation of human spliceosome associated factor 3, U4/U6 recycling protein
Symbol:	SART3
UniProt ID:	Q15020
Synonyms:	DSAP1; P100; p110; p110(nrb); RP11-13G14; TIP110
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)