

Product datasheet for **GA106548**

Human LPTM4A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LPTM4A CRISPRa kit - CRISPR gene activation of human lysosomal protein transmembrane 4 alpha
Symbol:	LPTM4A
Components:	GA106548G1, LPTM4A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106548G2, LPTM4A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106548G3, LPTM4A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q15012</u>
Synonyms:	HUMORF13; LPTM4; MBNT; Mtrp
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)