

Product datasheet for **GA106553**

Human **CLSTN3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	CLSTN3 CRISPRa kit - CRISPR gene activation of human calsynenin 3
Symbol:	CLSTN3
Components:	GA106553G1, CLSTN3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106553G2, CLSTN3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106553G3, CLSTN3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BQT9</u>
Synonyms:	alcbeta; CDHR14; CSTN3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)