


Product datasheet for **GA106554**

Human F115A (TCAF1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TCAF1 CRISPRa kit - CRISPR gene activation of human TRPM8 channel associated factor 1
Symbol:	TCAF1
Components:	GA106554G1, F115A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106554G2, F115A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106554G3, F115A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y4C2</u>
Synonyms:	FAM115A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)