

Product datasheet for **GA106557**

Human RIPOR2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	RIPOR2 CRISPRa kit - CRISPR gene activation of human RHO family interacting cell polarization regulator 2
Symbol:	RIPOR2
Components:	GA106557G1, RIPOR2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106557G2, RIPOR2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106557G3, RIPOR2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y4F9</u>
Synonyms:	C6orf32; DIFF40; DIFF48; MYONAP; PL48
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)