

Product datasheet for **GA106573**

Human PCLAF activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCLAF CRISPRa kit - CRISPR gene activation of human PCNA clamp associated factor
Symbol:	PCLAF
UniProt ID:	Q15004
Synonyms:	L5; NS5ATP9; OEATC; OEATC-1; OEATC1; p15(PAF); p15/PAF; p15PAF; PAF; PAF15
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)