

Product datasheet for **GA106576**

Human KIAA0195 (TMEM94) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM94 CRISPRa kit – CRISPR gene activation of human transmembrane protein 94
Symbol:	TMEM94
Components:	GA106576G1, KIAA0195 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106576G2, KIAA0195 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106576G3, KIAA0195 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q12767
Synonyms:	TMEM94
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

