

Product datasheet for **GA106578**

Human EIF4A3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EIF4A3 CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 4A3
Symbol:	EIF4A3
UniProt ID:	P38919
Synonyms:	DDX48; eIF4AIII; MUK34; NMP265; NUK34; RCPS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)