

Product datasheet for **GA106583**

Human FAM38A (PIEZO1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIEZO1 CRISPRa kit - CRISPR gene activation of human piezo type mechanosensitive ion channel component 1
Symbol:	PIEZO1
UniProt ID:	<u>Q92508</u>
Synonyms:	DHS; FAM38A; LMPH3; Mib
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)