

Product datasheet for **GA106598**

Human PHYHIP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PHYHIP CRISPRa kit - CRISPR gene activation of human phytanoyl-CoA 2-hydroxylase interacting protein
Symbol:	PHYHIP
Components:	GA106598G1, PHYHIP gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106598G2, PHYHIP gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106598G3, PHYHIP gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q92561</u>
Synonyms:	DYRK1AP3; PAHX-AP; PAHXAP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)