

Product datasheet for **GA106633**

Human KIAA0748 (TESPA1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TESPA1 CRISPRa kit - CRISPR gene activation of human thymocyte expressed, positive selection associated 1
Symbol:	TESPA1
Components:	GA106633G1, KIAA0748 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106633G2, KIAA0748 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106633G3, KIAA0748 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A2RU30
Synonyms:	HSPC257; KIAA0748
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)