

Product datasheet for **GA106650**

Human CEP170 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CEP170 CRISPRa kit - CRISPR gene activation of human centrosomal protein 170
Symbol:	CEP170
Components:	GA106650G1, CEP170 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106650G2, CEP170 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106650G3, CEP170 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5SW79</u>
Synonyms:	FAM68A; KAB; KIAA0470
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)