

Product datasheet for **GA106651**

Human LRIG2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRIG2 CRISPRa kit - CRISPR gene activation of human leucine rich repeats and immunoglobulin like domains 2
Symbol:	LRIG2
Components:	GA106651G1, LRIG2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA106651G2, LRIG2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA106651G3, LRIG2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q94898</u>
Synonyms:	LIG-2; LIG2; UFS2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)